

Chloroform Methanol Purification of Proteins

Nikolai Slavov

This protocol describes the extraction method used by Slavov *et al* (2014a,b). It is based on the classical method of Wessel and Flügge (1984) that has been used extensively for decades.

To 500 μ l protein sample (about 500 μ g protein) in a 15 ml falcon tube:

1. Add 2ml methanol and vortex thoroughly.
2. Add 500 μ l chloroform and vortex.
3. Add 1500 μ l water and vortex; the mixture become cloudy with precipitated protein flakes.
4. Centrifugation for 1 minute at $14,000 \times g$ resulted is three layers: a large aqueous layer on top, a circular flake of protein in the interphase, and a smaller chloroform layer at the bottom.
5. Remove top aqueous layer carefully, trying not to disturb the protein flake.
6. Add 2ml methanol and vortex.
7. Centrifuge the resulting mixture for 5 minutes at $20,000 \times g$ until the protein pellets.
8. Remove as much methanol as possible with care since the pellet is delicate.
9. Dry the protein pellet.

References

- Slavov N, Budnik B, Schwab D, Airoidi E, van Oudenaarden A (2014a) Constant Growth Rate Can Be Supported by Decreasing Energy Flux and Increasing Aerobic Glycolysis. *Cell Reports* **7**: 705 – 714
- Slavov N, Semrau S, Airoidi E, Budnik B, van Oudenaarden A (2014b) Variable stoichiometry among core ribosomal proteins. *arXiv preprint* **1**: arXiv:1406.0399
- Wessel D, Flügge UI (1984) A method for the quantitative recovery of protein in dilute solution in the presence of detergents and lipids. *Analytical biochemistry* **138**: 141–143